

# SAFEWAY

## Breath Alcohol Tester

### INSTRUCTIONS

#### Intended Use

The SAFEWAY is a pass/fail test to test for alcohol in human breath. The SAFEWAY is a disposable device designed for one-time use. The SAFEWAY alcohol test is a screening test that gives preliminary results.

#### Principle of Operation

The SAFEWAY uses a chemical reaction that changes color to indicate the presence of alcohol. When alcohol is present in the breath at a high enough concentration, the chemical reaction causes the yellow crystals to change to a green/blue or blue/green color. The SAFEWAY devices are calibrated to change color, indicating a positive test, at a specific breath alcohol level. The SAFEWAY comes in four different alcohol levels: .02, .04, .05, and .08%. The alcohol level for each device is printed on the label of the device.

#### Storage

Store at room temperature (50<sub>i</sub> to 85<sub>i</sub> F) away from direct sunlight.

#### Warnings, Precautions & Limitations

- Alcohol impairs judgment. Do not test yourself OR drive if you have been drinking. Do not use this test to determine if it is safe for you to drive. The actual result may be significantly higher or lower than indicated by this product. Someone who has not been drinking should help conduct the test and read the results.
- Keep out of the reach of children.
- Do not immerse in liquid.
- Do not inhale or eat the contents of the test. This product contains potassium dichromate, a hazardous chemical. If ingested, induce vomiting and contact your physician.
- The SAFEWAY alcohol tester is not intended to legally determine the presence of alcohol or concentration of alcohol in a person. This product should be used only as a screening device and is only an indication of the possible presence of alcohol in the blood of the user. The exact level of alcohol in the blood cannot be accurately determined by using this product. The actual result may be significantly higher or lower than indicated by this product.
- A positive result should be taken as a warning that a subject may have detectable alcohol in their system. A positive result should be confirmed by an evidentiary alcohol test before any legal or workplace actions are taken.
- There may be times that a person tests negative and later show that he or she is under the influence of alcohol or their judgment had been impaired by alcohol.
- For best results read test in indirect sunlight, incandescent or florescent lighting. The color change may be difficult to see under some streetlights and in dim light and poorly light areas.
- Do not use if glass tube is broken or if crystals are green or are not yellow.
- Do not interpret the test if color blind or visually impaired.
- Do not use after expiration date marked on package.

#### Disclaimer

The use of this device and the decisions of the user based on the results of the test are his/her sole responsibilities. The manufacturer, distributor or retailers are not responsible for the user's actions. The test is intended to give preliminary results only and does not determine intoxication levels. Positive test results suggest the presence of alcohol and should be confirmed with an evidential alcohol test.

Technical Assistance: Please call 866-216-8700 weekdays between 8 am and 5 pm for technical assistance, or email [info@safety-devices.com](mailto:info@safety-devices.com).

Made in the U.S.A.



Advanced Safety Devices, LLC.  
9140 Jordan avenue  
Chatworth, CA 91311

BR2500  
REV: 1  
2/1/12